

Unyvero

External Controls for LRT/Pneumonia Panel

unyvero



■ Simple and reliable external controls for your laboratory needs

The Unyvero LRT/Pneumonia Control Kit includes 12 single-use vials of ready-to-use qualitative external controls with coverage for all targets on the Unyvero LRT Cartridge. Each kit contains:

- 4 tubes x Unyvero LRT/Pneumonia Positive A
- 4 tubes x Unyvero LRT/Pneumonia Positive B
- 4 tubes x Unyvero LRT/Pneumonia Negative

Ordering information: Contact MMQCI for all orders of Unyvero LRT/Pneumonia Control Panel M404. Submit orders by phone, fax or email:

phone: **1-207-885-1072 ext. 1**

fax: 1-207-885-1079

email: **info@mmqci.com**

Storage Conditions: Refrigerated 2 – 8°C. Unopened tubes are stable through the expiration date printed on the kit label when appropriately stored.

MMQCI Part Number: M404

Price: \$427 per kit

Questions: Contact MMQCI at **1-207-885-1072**



Curetis USA
619-452-3644
info_us@curetis.com
www.curetis.com